

Group Primer Pairs Element

Select groups of primer pairs, which can be simultaneously used in one reaction tube.

The primers must be supplied in the following order: pair1_direct_primer, pair1_reverse_primer, pair2_direct_primer, pair2_reverse_primer, etc.

Parameters in GUI

Parameter	Description	Default value
Output report file	Path to the report output file.	

Parameters in Workflow File

Type: primers-grouper

Parameter	Parameter in the GUI	Type
output-file	Output report file	string

Input/Output Ports

The element has 1 *input port*:

Name in GUI: *Primer pairs*

Name in Workflow File: in-sequence

Slots:

Slot In GUI	Slot in Workflow File	Type
Sequence	sequence	<i>sequence</i>