


# Search in Sequence

To search for a pattern(s) in a sequence go to the *Search in Sequence* tab of the the *Options Panel* in the *Sequence View*.

Input the value you want to search in the text field and click the *Search* button. To search multiple patterns input the patterns separated by a new line in the pattern text field. To add a new line symbol *Ctrl+Enter* may be used. You can input the value as sequence or name of the sequence in the FASTA format and sequence after that.



By default, **misc\_feature** annotations are created for regions that exactly match the pattern. Find below the description of the available settings.

- [Load Patterns from File](#)
- [Search Algorithm](#)
- [Search in](#)
- [Other Settings](#)
- [Annotations Settings](#)