


# General



The following settings are available on the tab:

*Language of User Interface (applied after restart)* — here you can select UGENE localization. Currently, available localizations are EN and RU. The default value (*Autodetection*) specifies that UGENE should use the operating system regional options to select the localization. This setting is applied only after UGENE is reopened.

*Appearance* — defines the appearance of the application.

*Open last project at startup* — if the option is checked, the last project is opened when UGENE is started. Also, you can choose default settings for a saving project.

*Enable statistical reports collecting* — collects information about UGENE usage and sends it to the UGENE team to help improve the application.



The collected information includes:

1. System info: UGENE version, OS name, Qt version, etc.
2. Counters info: number of launches of certain tasks (e.g. HMM search, MUSCLE align).

The collected information DOESN'T include any personal data.

*Check for updates when UGENE starts* - checks for updates when UGENE starts.

*Default settings* — this option resets the default settings on the next run.