

Workflow Elements

This section contains detailed description of all workflow elements presented in the Workflow Designer.

For each element you can find:

- Description of the parameters used in the GUI
- Corresponding parameters names used in a workflow file
- Information about input and output ports

The type of a parameter can be one of the following:

string

A string.

numeric

A number.

boolean

A boolean data type. Available values are: true / false, 0 / 1 and yes / no.

A port's slot type can be one of the following:

sequence

Biological sequence

msa

Multiple sequence alignment

text

A text

annotation-table

Table of annotations

annotation-table-list

A list of different tables of annotations

ebwt-index

Bowtie index

hmm2-profile

A HMM profile of HMMER2 package

fmatrix

Frequency matrix

wmatrix

Weight matrix

sitecon-model

SITECON model


assembly

Assembly

variation

Variation track

To search an element use the name filter or press the *Ctrl+F* shortcut that moves you to the name filter also:



- Data Readers
 - File List Element
 - Read Alignment Element
 - Read Annotations Element
 - Read Assembly Element
 - Read from DAS Element
 - Read from Remote Database Element
 - Read Plain Text Element
 - Read Sequence Element
 - Read Variations Element
- Data Writers
 - Write Alignment Element
 - Write Annotations Element
 - Write Assembly Element
 - Write FASTA Element
 - Write Plain Text Element
 - Write Sequence Element
 - Write Variations Element
- Data Flow
 - Filter Element
 - Grouper Element
 - Multiplexer Element
 - Sequence Marker Element
- Basic Analysis
 - Amino Translations Element
 - Annotate with DAS Element
 - Annotate with UQL Element
 - CD-Search Element
 - Collocation Search Element
 - Export PHRED Qualities Element
 - Fetch Sequences by ID From Annotation Element
 - Filter Annotation by Name Element
 - Filter Annotations by Qualifier
 - Find Pattern Element
 - Find Repeats Element
 - Gene-by-gene approach report
 - Get Sequences by Annotations Element
 - Import PHRED Qualities Element
 - Local BLAST Search Element
 - Local BLAST+ Search Element
 - Merge Annotations Element
 - ORF Marker Element
 - Remote BLAST Element
 - Remove Duplicates in BAM Files Element
 - Smith-Waterman Search Element
- Data Converters
 - Convert bedGraph Files to bigWig Element
 - Convert Text to Sequence Element
 - File Format Conversion Element
 - Reverse Complement Element
 - Split Assembly into Sequences Element
- DNA Assembly
 - Align reads with BWA-MEM
 - Assembly Sequences with CAP3
 - Extract Consensus from Assembly
- HMMER2 Tools
 - HMM Build Element
 - HMM Search Element
 - Read HMM Profile Element
 - Write HMM Profile Element
- HMMER3 Tools
 - HMM3 Build Element
 - HMM3 Search Element
 - Read HMM3 Profile
 - Write HMM3 Profile
- Multiple Sequence Alignment
 - Align Profile to Profile with MUSCLE Element
 - Align with ClustalO Element
 - Align with ClustalW Element
 - Align with Kalign Element
 - Align with MAFFT Element
 - Align with MUSCLE Element
 - Align with T-Coffee Element
 - Extract Consensus from Alignment
 - Join Sequences into Alignment Element
 - Split Alignment into Sequences Element
- NGS: Basic
 - CASAVA FASTQ Filter Element

- FASTQ Quality Trimmer Element
- Filter BAM/SAM Files Element
- Genome Coverage Element
- Merge BAM Files Element
- Slopbed Element
- Sort BAM Files Element
- NGS: ChIP-Seq Analysis
 - Annotate Peaks with peak2gene Element
 - Build Conservation Plot Element
 - Collect Motifs with SeqPos Element
 - Conduct GO Element
 - Create CEAS Report Element
 - Find Peaks with MACS Element
- NGS: RNA-Seq Analysis
 - Assembly Transcripts with Cufflinks Element
 - Extract Transcript Sequences with gffread Element
 - Find Splice Junction with TopHat Element
 - Merge Assemblies with Cuffmerge Element
 - Test for Diff. Expression with Cuffdiff Element
- NGS: Variant Calling
 - Call Variants with SAMtools Element
 - Create VCF consensus
- SNP Annotation
 - Annotate variations with SNPToolbox Element
 - Detect Transcription Factors with rSNP-Tools Element
 - Determine SNP effect on TATA-boxes Element
 - ProtStability1D Element
 - ProtStability3D Element
 - SNP Chip Tools Element
 - SNP Effect on PDB sites Element
 - Write SNP Report Element
- Transcription Factor
 - Build Frequency Matrix Element
 - Build SITECON Model Element
 - Build Weight Matrix Element
 - Convert Frequency Matrix Element
 - Read Frequency Matrix Element
 - Read SITECON Model Element
 - Read Weight Matrix Element
 - Search for TFBS with SITECON Element
 - Search for TFBS with Weight Matrix Element
 - Write Frequency Matrix Element
 - Write SITECON Model Element
 - Write Weight Matrix Element
- Utils
 - DNA Statistics Element
 - Generate DNA Element
- Custom Elements With Script
 - CASAVA FASTQ Filter Script Element
 - Dump Sequence Info Element
 - FASTQ Trimmer Element
 - LinkData Fetch Element
 - Quality Filter Element