


Managing Zoom View Rows

The Zoom View contains not more than 20 rows by default. The rest rows are available by scrolling.

To change this behavior use the *Manage Rows in Zoom View* menu button on a sequence toolbar:



When the *Show All Rows* item is checked all available annotations are always shown. You can also add rows by selecting the *+5 Rows* and *+1 Row* items and remove rows by selecting the *-5 Rows* and *-1 Row* items. To restore the default number of rows select the *Reset Rows Number* item.

See also:

- [Navigating Sequence zoom view using Sequence overview](#)
- [Zooming Sequence](#)
- [Creating New Ruler](#)
- [Manipulating Annotations](#)