

Creating Document

To create a new sequence file from a text, select the *File New document from text* main menu item or *Create Sequence* on the main screen.

The *Create Sequence* dialog appears:

Create Sequence

Paste data here

ACGT

☒ **Custom settings**

Alphabet: Standard DNA ▼

☒ Skip unknown symbols

☐ Replace unknown symbols with

Save sequence to file ...

File format FASTA ▼

Sequence name Sequence

☐ Save file immediately

Help Cancel Create

You can input the created sequence to the *Paste data here* field. You can type or paste sequences in FASTA or text format.

The following *Custom settings* are available:

Alphabet — here you can select the alphabet:

☒ **Custom settings**

Alphabet: Standard DNA ▼

- Standard DNA
- Standard RNA
- Extended DNA
- Extended RNA
- Standard amino acid
- Extended amino acid
- Raw

Save sequence to file ...

File format FASTA ▼

Sequence name Sequence

The following alphabets are available: [Standard DNA](#), [Standard RNA](#), [Extended DNA](#), [Extended RNA](#), [Standard amino](#), [Extended amino](#), Raw.

Skip unknown symbols / Replace unknown symbols with — you can select either to skip unknown input symbols or to replace them with the specified symbol.

Save sequence to file — a location of the created document.

File format — a format of the created document. Currently available formats are FASTA and Genbank.

Sequence name — a name of the sequence in the created document.

Save file immediately — check this option if you want to save the document immediately after the *Create* button is pressed.

The created document will be added to the current project and opened in the *Sequence View*.