

# Captions on Annotations

It is possible to show a value of a qualifier of an annotation instead of the annotation type name in the *Sequence Zoom View*. To enable this option for an annotation type check the *Show value of qualifier* check box and input the values of the required qualifiers in the text field nearby this check box. See the image below.

The screenshot displays the NC\_001363 [d] sequence viewer. The main window shows a genomic track with a scale from 1 to 5k. A specific region is highlighted with a blue arrow indicating a 589 bp segment. Below the scale, the DNA sequence is shown with various annotations. A red arrow points to a specific annotation in the track.

The 'Annotations Highlighting' panel on the right allows users to select an annotation name and configure its display. The 'Show value of qualifier' checkbox is checked, and the 'note' text field is visible. A red arrow points to this checkbox.

| Annotation   | Color       |
|--------------|-------------|
| BadI         | Light Green |
| BaeGI        | Orange      |
| BaeI         | Purple      |
| CDS          | Green       |
| SacII        | Light Blue  |
| comment      | Pink        |
| misc_feature | Brown       |
| source       | Grey        |

| Name                      | Type          | Value                   |
|---------------------------|---------------|-------------------------|
| NC_001363 features [m...] |               |                         |
| source (0, 1)             |               |                         |
| misc_feature (0, 2)       |               |                         |
| misc_feature              | Misc. Feature | 5245..5833              |
| note                      |               | Show value of qualifier |
| misc_feature              | Misc. Feature | 2..590                  |
| comment (0, 1)            |               |                         |
| CDS (0, 4)                |               |                         |
| CDS                       | CDS           | 5048..5203              |
| CDS                       | CDS           | 3875..4999              |

If you input several qualifiers names (separated by comma), then the first found qualifier is taken into account and shown on the annotation.